

# Analysis of Directly Cast Metallurgical Slags

## ARL OPTIM'X WDXRF sequential configuration

### Key Words

- ARL OPTIM'X
- Iron & Steel
- Slags
- WDXRF
- X-Ray Fluorescence

### Introduction

Slags originate from various stages in the iron & steel process, e.g. blast furnace, converter, basic oxygen furnace (BOF also referred to as BOS = basic oxygen steelmaking), electric arc furnace or ladle.

In the blast furnace, slag is formed from impurities in the iron ores (known as the gangue), the flux and coke ashes. It is a complex silicate of aluminum, calcium and magnesium containing small quantities of oxides of manganese and iron as well as calcium sulfide. Slag has a double role: it permits removal of the gangue thanks to its fusibility and fluidity, it also allows exchange reaction with the liquid metal and permits a control of the process in order for the desirable elements to stay in the melt while the others are removed. As an example, in an electric arc furnace the slag formation process can be controlled by adding oxygen, carbon and slag formers to the melt. This will promote formation of CO instead of MnO and FeO and help keep these elements in their metallic form in the melt. The basic slag formers like lime (CaO) and magnesia (MgO) will help neutralize the acidity of the slag in order to save the refractory bricks of the furnace.

### Instrumentation

The Thermo Scientific ARL OPTIM'X is a WDXRF instrument designed for ease of use, minimal operating and maintenance costs. An ARL OPTIM'X XRF spectrometer configured as a sequential unit was used to obtain the results presented in this application note. The instrument was fitted with a SmartGonio™ covering elements from F



(Z=9) to U (Z=92). The ARL OPTIM'X XRF spectrometers are fitted with a low power Rh anode X-ray tube and the geometry of the instrument is optimized to provide the highest sensitivity. The instrument does not require external or internal water cooling, and has superior resolution, precision, short and long term

stability. It can analyse Na and Mg without problem or even F if necessary. Ease of operation is obtained through OXSAS, our modern, powerful and user-friendly software supporting instrument operation and data handling.

### Sample preparation

While still hot the slag sample is directly cast into a steel ring using a ceramic mold specially designed for this



Figure 1: Directly cast slags

purpose. This direct cast sample can be analyzed in the XRF instrument without further preparation (Figure 1).

### Calibration and results

Ten secondary slag standards were used for calibration of the ARL OPTIM'X. These standard samples allow the concentration ranges shown in Table 1 to be covered. A working curve is established for each element using the Multi-Variable-Regression incorporated in the OXSAS software package. The Standard Error of Estimate (SEE) is a measure of the accuracy of analysis. It is the average error between the certified concentration of the standard samples and the calibration curve of a given oxide.

ELEMENTS/ OXIDES	ANALYTICAL DEVICE	CALIBRATION RANGES	SEE (%)
Al <sub>2</sub> O <sub>3</sub>	SmartGonio™	5 % - 12 %	0.074 %
MgO	SmartGonio™	2.4 % - 9 %	0.14 %
SiO <sub>2</sub>	SmartGonio™	34 % - 38 %	0.23 %
CaO	SmartGonio™	32 % - 47 %	0.24 %
Cr <sub>2</sub> O <sub>3</sub>	SmartGonio™	1.6 % - 7.8 %	0.16 %

Table 1: Summary of performance

SEE = Standard error of estimate : it is a measure of the accuracy

## Stability tests

In order to show the excellent repeatability of the ARL OPTIM'X for the analysis of slags both in short and long term, stability tests were performed. For short term repeatability 10 consecutive measurements were performed on one of the cast slag samples. The repeatability test was done twice using counting times of 20s and 40s per element. Average concentration and standard deviations are shown in Table 2.

ELEMENT/OXIDE	ANALYTICAL DEVICE	AVERAGE CONCENTRATION %	COUNTING TIME	STD.DEV. %	COUNTING TIME	STD.DEV. %
Al <sub>2</sub> O <sub>3</sub>	SmartGonio™	9.74	20 s	0.034	40 s	0.027
CaO	SmartGonio™	44.1	20 s	0.034	40 s	0.033
Cr <sub>2</sub> O <sub>3</sub>	SmartGonio™	4.59	20 s	0.018	40 s	0.017
MgO	SmartGonio™	3.17	20 s	0.018	40 s	0.013
SiO <sub>2</sub>	SmartGonio™	36.2	20 s	0.041	40 s	0.026
Total counting time		100 s	200 s			

Table 2: Results of a repeatability tests (10 runs) using counting times of 20s and 40s per element

## Conclusion

Various elements in directly cast slags can be successfully analyzed using the ARL OPTIM'X low power WDXRF instrument. The results obtained on cast slags are very similar to those obtained on slags prepared through grinding and pressing (see application notes nr. 41703 and 41704)

The limitation in term of accuracy is not due to the instrument but to the sample preparation and the variety of slags used in this test. In order to get the best accuracy with slags of different origins it is necessary to prepare the samples as fused beads with Lithium tetraborate (see application note nr. 41702) in order to avoid grain size effects and mineralogical effects.

Good repeatability and reproducibility is obtained with 100s of total counting time with the SmartGonio™. If better results are required for any element the counting time for that particular element can be increased. Additionally two fixed channels can be fitted alongside the SmartGonio™ of the ARL OPTIM'X in order to improve the performance for two specific elements when required.

To know more on our complete X-ray product portfolio, please visit [www.thermo.com/xray](http://www.thermo.com/xray).

## Acknowledgments

We would like to thank Mr. Khalid El Mansour and Mr. Jean-Claude Lemoine of Arcelor for providing the cast slag samples and their concentrations originating from the Ardoise plant.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

### Africa-Other

+27 11 570 1840

### Australia

+61 2 8844 9500

### Austria

+43 1 333 50 34 0

### Belgium

+32 2 482 30 30

### Canada

+1 800 530 8447

### China

+86 10 8419 3588

### Denmark

+45 70 23 62 60

### Europe-Other

+43 1 333 50 34 0

### Finland/Norway/ Sweden

+46 8 556 468 00

### France

+33 1 60 92 48 00

### Germany

+49 6103 408 1014

### India

+91 22 6742 9434

### Italy

+39 02 950 591

### Japan

+81 45 453 9100

### Latin America

+1 608 276 5659

### Middle East

+43 1 333 50 34 0

### Netherlands

+31 76 579 55 55

### South Africa

+27 11 570 1840

### Spain

+34 914 845 965

### Switzerland

+41 21 694 71 11

### UK

+44 1442 233555

### USA

+1 800 532 4752

[www.thermo.com](http://www.thermo.com)



Thermo Fisher Scientific (Ecuublens) SARL, Switzerland is ISO certified.

AN41707\_E 08/09C

©2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.