

Impurity Profiling of Carbamazepine by HPLC/UV

Terry Zhang, Guifeng Jiang, Thermo Fisher Scientific, San Jose, CA, USA

Key Words

- Accela
- Hypersil GOLD
- Carbamazepine
- Drug Analysis

Goal

Develop a simple, accurate and robust HPLC/UV method for quantitative impurity analysis of carbamazepine.

Introduction

Carbamazepine is a medication indicated for use in epilepsy, trigeminal neuralgia and bipolar disorder.¹ This widely prescribed tricyclic anticonvulsant is administered orally in usual doses of 400–1200 mg/day and is available in several commercial tablet, capsule and suspension forms.¹ Due to its narrow therapeutic range, tight control and routine monitoring of plasma carbamazepine concentrations are required to optimize treatment while preventing adverse drug effects.¹ Safe and effective therapy necessitates reliable and accurate methods for assessing all aspects of drug quality, including purity, strength and stability.

Synthesis-related organic impurities are usually present in bulk forms as well as in pharmaceutical formulations of carbamazepine. Pharmacopoeias set strict standards for the purity of carbamazepine to ensure drug efficacy and safety. The United States Pharmacopeia (USP) establishes 0.2% as the maximum limit for any individual impurity and 0.5% as the total of all impurities relative to the active pharmaceutical ingredient.² The International Conference on Harmonization (ICH) guidelines recommend the reporting of any impurity above a level of 0.05% and the identification of all impurities above a level of 0.10% in new drug substances that are administered at <2 g/day doses.³ The USP describes a quantitative HPLC/UV procedure for impurity profiling of carbamazepine that utilizes a 4.6 mm × 250 mm column with L10 packing (5–10 μm silica particles with cyano bonded phases) and a mobile phase consisting of water, methanol, tetrahydrofuran, formic acid and triethylamine. However, the utility of this method for routine quality control analysis is limited by the complexity of the mobile phase, which could lead to poor reproducibility, as well as by the difficulty in achieving the required resolution ($R_s \geq 1.7$) between carbamazepine and the impurity 10,11-dihydroxycarbamazepine (carbamazepine related compound A), which have comparable polarities and differ widely in their concentrations. Routine drug purity analysis requires simple and robust analytical techniques that deliver exceptional resolution, sensitivity, and accuracy.



The Thermo Scientific Accela UHPLC system and Thermo Scientific Hypersil GOLD columns together enable sensitive, accurate and robust LC analysis for pharmacopeia methods and maximize flexibility for method development. The Accela™ 1250 Pump delivers precise flows and accurate gradients at an expansive range of flow rates (up to 2 mL/min) and pressures (up to 1250 bar), and the Accela PDA detector offers up to five times the sensitivity of conventional photodiode array detectors. Hypersil GOLD™ columns enhance peak shapes, resolution and reproducibility through high purity silica technology and robust bonding and end-capping procedures. In this note, we describe a simple HPLC/UV assay for accurate and robust impurity analysis of carbamazepine using the Accela UHPLC system and a high performance Hypersil GOLD CN (cyano) column, and demonstrate that the method meets USP requirements and ICH guidelines.



Materials and Methods

Sample Preparation

Carbamazepine and related impurities A, B, D and F were purchased from Sigma Aldrich. Stock solutions of 2 mg/mL were prepared in methanol. Solutions of 100 µg/mL concentration were prepared by diluting in methanol and were used for method development. Calibration solutions were prepared by serial dilution of the stock solutions in methanol, at concentrations of 250 ng/mL–400 µg/mL for each impurity and 5 µg/mL–1.5 mg/mL for carbamazepine. A mixture containing 1 mg/mL of carbamazepine and 500 ng/mL of each impurity was prepared to examine method suitability at a 0.05% reporting threshold. For assessments of system suitability and quantification, a mixture of 1.5 mg/mL of carbamazepine with 3 µg/mL of each impurity was used.

HPLC/UV Analysis

Instrumentation

HPLC separations were performed on an Accela 1250 UHPLC system and an Accela autosampler (Thermo Fisher Scientific, San Jose, CA, USA). UV absorbance was monitored at 211 nm at 20 Hz using an Accela PDA detector (Thermo Fisher Scientific, San Jose, CA, USA).

LC Parameters:

Columns:	Thermo Fisher Scientific Hypersil GOLD aQ column (2.1 × 100 mm, 1.9 µm particle size) and Hypersil GOLD CN column (2.1 × 150 mm, 3.0 µm particle size)			
Mobile Phase:	A: Water B: Acetonitrile			
Column Temperature:	30 °C			
Sample Injection Volume:	1 µL			
Needle Wash:	80:20 (v/v) acetonitrile:water			
Gradient:	Time	A %	B %	µL/min
	0.00	88.0	12.0	400.0
	10.00	88.0	12.0	400.0
	14.00	75.0	25.0	400.0
	28.00	70.0	30.0	400.0
	19.00	15.0	85.0	400.0
	30.00	15.0	85.0	400.0
	31.00	88.0	12.0	400.0

Results and Discussion

HPLC Separation of Carbamazepine and Related Impurities

The major organic impurities that arise from the synthesis of carbamazepine are 10,11-dihydrocarbamazepine (impurity A, also referred to as USP carbamazepine related compound A), 9-methylacridine (impurity B), *N*-carbamoylcarbamazepine (impurity C), iminostilbene (impurity D), iminodibenzyl (impurity E), and 5-chlorocarbonyliminostilbene (impurity F). An HPLC/UV assay that can accurately and reproducibly detect, identify and quantitate these compounds at trace levels (as low as 0.05%–0.20% relative to the parent drug carbamazepine) is required to comply with USP requirements and ICH guidelines.

The USP chromatographic method for determining carbamazepine purity specifies the use of a 4.6 mm × 250 mm L10 column (cyano column with 5–10 µm silica particles), a mobile phase composed of 1000 mL water/methanol/tetrahydrofuran (85:12:3) with 0.22 mL formic acid and 0.5 mL triethylamine, a flow rate of approximately 1.5 mL/min, and detection at 230 nm. Drawbacks of this assay for QA/QC applications include the use of a complex mobile phase with an ion-pairing agent, which could lead to operational difficulties and poor reproducibility, and the difficulty in meeting the resolution requirement ($R_s \geq 1.7$) for carbamazepine and impurity A. Resolution and method performance may be improved by using high efficiency HPLC columns with smaller (≤ 3 µm) particles and by modifying the mobile phase.

Using the Accela 1250 UHPLC system, a Hypersil GOLD aQ polar endcapped C18 column (1.9 µm, 2.1 × 100 mm) and a simple acetonitrile/water gradient, a standard mixture containing carbamazepine and impurities A, B, D and F at 100 µg/mL concentrations was separated and detected within 23 minutes (Figure 1a). Detection at 211 nm was determined to provide the optimum response for this analysis; acetonitrile exhibits little absorption at this wavelength. All analytes were baseline-resolved and eluted in the order of carbamazepine, impurity A, impurity D, impurity B and impurity F. Symmetric chromatographic peaks were obtained and a resolution of 2.02 between carbamazepine and impurity A was achieved. However, resolution between this critical pair was considerably reduced when the concentration of impurity A was 0.1% of the parent drug (Figure 1b). Enhanced resolution is required to ensure compliance with USP requirements.




Figure 1a: HPLC separation of a mixture of carbamazepine and its related impurities (100 µg/mL) using a 1.9 µm Hypersil GOLD aQ column, 2.1 mm × 100 mm. Mobile phase: A – Water; B – Acetonitrile. UV detection: 211 nm.

Figure 1b: Impurity A at a level of 0.1% is partially masked by the parent drug. Higher resolution (>1.7) between this critical pair is required to meet USP requirements.

The Hypersil GOLD CN (cyano) column offers alternative selectivity in reversed phase chromatography with lower hydrophobicity compared to C18 alkyl chain phases. The optimal wavelength of 211 nm was determined based on the best signal-to-noise ratio (S/N) for all analytes (Figure 2). Baseline separation of carbamazepine and the impurities at 100 µg/mL concentrations was achieved in 22 minutes using a Hypersil CN column (3 µm, 2.1 × 150 mm) and a simple acetonitrile/water gradient (Figure 3). Carbamazepine, impurity A and impurity D were less retained on the cyano phase and the elution order of carbamazepine and impurity A was reversed with this column. Excellent UV responses of similar magnitude were obtained for all analytes (Figure 3). The impurity responses were calculated relative to carbamazepine and were applied later for impurity analysis and quantification of a standard carbamazepine solution. The inset table of Figure 3 shows the responses of the impurities relative to Carbamazepine. Maximizing peak heights leads to sensitivity gains and is especially beneficial in trace analysis. Compared to Hypersil aQ column, the Hypersil CN column notably improved the peak shapes of carbamazepine and impurity A to enable enhanced resolution and sensitivity at ICH reporting and identification thresholds (0.05% and 0.10%, respectively).




Figure 2: Selection of the wavelength for quantitation was based on the best S/N for all analytes




Figure 3: HPLC separation of a mixture of carbamazepine and its related impurities (100 µg/mL) using a 3 µm Hypersil GOLD CN column, 2.1 mm × 150 mm. Mobile phase: A – Water; B – Acetonitrile. UV detection: 211 nm.

A standard solution of 1 mg/mL carbamazepine was analyzed and 600 ng/mL of impurity A (0.06%) was detected (Figure 4). The resolution between the carbamazepine and impurity A peaks was 2.28, which exceeds the minimum resolution specified by USP (Figure 4). Figure 5 shows the chromatogram of the standard 1 mg/mL carbamazepine solution spiked with 0.05% of impurities A, B, D and F. All impurities were easily detected at these trace amounts, and the stronger UV response exhibited by impurity A was attributed to the 0.06% that was already present in the carbamazepine standard. As the majority of aqueous and matrix contaminants usually elute early at void volume, elution of the first analyte of interest after the sixth minute ensures a robust quantitation method.

Method Validation

System suitability was investigated by analyzing six replicate injections of a solution containing 1500 µg/mL of carbamazepine and 3 µg/mL of each impurity, corresponding to 0.2% of the parent drug. Retention time RSDs ranged from 0.01–0.02%, while peak area RSDs ranged from 1–2% (Table 1), indicating excellent method reproducibility, particularly of the LC pump. According to Equation 1, signal-to-noise ratios of ≥ 25 are required to comply with the $\leq 2\%$ RSD criterion specified by the USP method:

$$\% \text{ RSD} \approx 50/(S/N) \quad (\text{Equation 1})$$




Figure 4: 1.0 mg/mL of carbamazepine standard was analyzed and 0.06% impurity A (600 ng/mL) was detected. Excellent resolution is achieved between carbamazepine and impurity A on the Hypersil GOLD CN column. The resolution between this critical pair is 2.28, which exceeds the USP criterion of >1.7 .




Figure 5: Impurities are baseline-separated and easily detected at the ICH reporting threshold (0.05%). Inset is the zoomed chromatogram. Mixture contains 1.0 mg/mL of carbamazepine and 500 ng/mL each of the impurities A, B, D, and F. Column: 3 µm Hypersil GOLD CN, 2.1mm × 150 mm. Mobile phase: A – Water; B – Acetonitrile. UV detection: 211 nm.

Compound	Concentration µg/mL	n	RT RSD%	Area RSD%	Mean Rs	Mean S/N	CAA*
Impurity_A	3	6	0.01	2	NA	73.4	111
Carbamazepine	1500	6	0.01	1	2.2	19720.8	NA
Impurity_D	3	6	0.01	2	9.1	44.1	108
Impurity_B	3	6	0.02	2	11.0	28.5	113
Impurity_F	3	6	0.01	1	2.2	34.6	102

Table 1: Method validation parameters were assessed by analyzing six replicate injections of each compound. Column: 3 µm Hypersil GOLD CN, 2.1 mm × 150 mm. Mobile phase: A – Water; B – Acetonitrile. UV detection: 211 nm.

CAA*: corrected area % accuracy

Mean signal-to-noise ratios of the impurities at 0.2% levels relative to the parent drug ranged from 28.5 to 73.4 (Table 1), ensuring method reliability for quantification and compliance with USP requirements. Mean resolutions ranged from 2.2 to 11.0, and the mean resolution of 2.2 between carbamazepine and impurity A exceeded the USP resolution requirement of 1.7 (Table 1).

Quantitative accuracy for all impurities was excellent, ranging from 102 to 111% (Table 1). Table 2 shows that excellent linearity in detector response was observed over the range of 0.25–400 µg/mL (ppm) for all impurities and 5–1500 µg/mL (ppm) for carbamazepine, with correlation coefficients ≥0.997 for all analytes.

Conclusion

A simple, accurate and robust quantitative HPLC/UV assay for carbamazepine purity assessment was developed. Separation and detection of carbamazepine and four of its related synthetic impurities was achieved within 22 minutes using the Accela 1250 UHPLC system, a Hypersil GOLD CN column and a simple water/acetonitrile gradient. The highly efficient and selective Hypersil GOLD CN column enhanced resolution between carbamazepine and impurity A, provided sensitivity gains through improved peak shapes, and eliminated the need for ion-pairing agents. Trace levels of impurities (0.05% of parent drug) were easily detected. Impurities at levels of 0.2% of carbamazepine were quantifiable with ≤2% RSD. This simple HPLC/UV assay meets USP criteria and ICH guidelines and provides a powerful alternative to the HPLC-based procedures currently recommended by regulatory agencies for QA/QC of carbamazepine purity.

www.thermoscientific.com

Legal Notices: ©2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Compound	Concentration µg/mL	Correlation Coefficients
Impurity_A	0.25–400	1.0000
Carbamazepine	5–1500	0.9973
Impurity_D	0.25–400	0.9995
Impurity_B	0.25–400	0.9999
Impurity_F	0.25–400	0.9998

Table 2: Excellent linearity of detector response was achieved, with correlation coefficients ≥0.997 for all compounds

References

- <http://www.merckmanuals.com/professional/lexicomp/carbamazepine.html>
- U.S. Pharmacopeia Monograph for Carbamazepine, USP32-NF27, page 1784.
- International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use, Impurities in New Drug Substances Q3A(R2), 2006.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other
+27 11 570 1840

Australia
+61 3 9757 4300

Austria
+43 1 333 50 34 0

Belgium
+32 53 73 42 41

Canada
+1 800 530 8447

China
+86 10 8419 3588

Denmark
+45 70 23 62 60

Europe-Other
+43 1 333 50 34 0

Finland/Norway/Sweden
+46 8 556 468 00

France
+33 1 60 92 48 00

Germany
+49 6103 408 1014

India
+91 22 6742 9434

Italy
+39 02 950 591

Japan
+81 45 453 9100

Latin America
+1 561 688 8700

Middle East
+43 1 333 50 34 0

Netherlands
+31 76 579 55 55

New Zealand
+64 9 980 6700

Russia/CIS
+43 1 333 50 34 0

South Africa
+27 11 570 1840

Spain
+34 914 845 965

Switzerland
+41 61 716 77 00

UK
+44 1442 233555

USA
+1 800 532 4752



Thermo Fisher Scientific, San Jose, CA USA is ISO Certified.

AN52049_E 03/11M