



Monitoring Protein PEGylation with Ion Exchange Chromatography

Peter Yu, Deanna Hurum, Leo Wang, Terry Zhang,
and Jeffrey Rohrer

245th ACS National Meeting & Exposition


April 9, 2013

OT70618_E 04/13S

PEGylation of Therapeutic Proteins

- Polyethylene glycol (PEG)ylation is used to modify therapeutic proteins.
 - Improved stability
 - Reduced circulatory clearance
 - Potential for smaller doses/greater time between doses
- It is used in numerous pharmaceutical agents for delivery.
- Most often, *N*-Hydroxysuccinimide (NHS) is used to react with lysine amines to form a covalent attachment with the protein.

Potential Challenges




- Positional Isomerism
- Site Specification
- Separation of Complex Mixtures

Ion Exchange and PEGylated proteins

- Two model proteins:
 - Lysozyme
 - RNase A
- A small but commercially available PEG was used to modify these model proteins.
- The reaction conditions, varying both time and the reactant ratio, were investigated.

Reagents

- Branched Amine-Reactive PEG



- Lysozyme—six potential lysine sites, 14.3 kD
- RNase A—ten potential lysine sites, 13.7 kD

Sample Preparation

- Dissolve 1 mg protein in 1 mL of phosphate-buffered saline.
- Dispense appropriate volume of PEGylation reagent into protein solution to achieve the desired reactant molar ratios.
- Incubate at room temperature for 30 min.
- Desalt using 0.5 mL Thermo Scientific™ Zeba™ Spin Desalting Columns to remove buffers and unreacted PEGylation reagent.

Instrumentation

- Thermo Scientific™ Dionex™ UltiMate™ 3000 BioRS System consisting of:
 - SRD-3600 Integrated Solvent and Degasser Rack
 - DGP-3600RS Biocompatible Dual-Gradient Rapid Separation Pump
 - TCC-3000SD Thermostatted Column Compartment
 - WPS-3000TBRS Thermostatted Biocompatible Rapid Separation Autosampler
 - VWD-3400 Variable Wavelength Detector


Chromatographic Method

- Column: Thermo Scientific™ ProPac™ SCX-10G Guard (2 × 50 mm), ProPac SCX-10 Analytical (2 × 250 mm)
- Mobile Phases: A) 20 mM MES, pH 6.1
B) 20 mM MES, 1M NaCl, pH 6.1
- Gradient: 0–30% B from 0 to 30 min, equilibration at 0% B for 10 min before injection
- Flow Rate: 0.25 mL/min
- Inj. Volume: 10 µL
- Detection: UV, 280 nm

Mass Spectrometry



- The PEGylation reaction mixture was desalted and injected on a Thermo Scientific™ BioBasic™ 8 HPLC Column (5 μ m, 100 \times 1.0 mm).
- MS confirmation was performed on a Thermo Scientific™ Q Exactive™ Orbitrap™ mass spectrometer with 140,000 mass resolution. The Q Exactive was operated in positive full-scan mode with scan range from 1000 to 3000 m/z.
- Mass spectra deconvolution was accomplished using the Xtract Algorithm of the Thermo Scientific™ Xcalibur™ software.

Literature Separation of PEGylated Lysozyme




Biotechnol. J. **2010**, 5, 477–483.

PEGylated Lysozyme Separation from a Vendor's Application Note



ProPac IEX Resin with Grafted Polymeric Functionalized Chains



Changes in Lysozyme PEGylation with Reaction Time

Columns: ProPac SCX-10G Guard, 2 × 50 mm
ProPac SCX-10 Analytical, 2 × 250 mm


Mobile Phases: A) 20 mM MES, pH 6.1
B) 20 mM MES, 1M NaCl, pH 6.1

Gradient: 0–30% B from 0 to 30 min
Equilibration at 0% B for 10 min before injection


Flow Rate: 0.25 mL/min

Inj. Volume: 10 µL


Samples: Lysozyme (black)
Lysozyme reacted with PEG for 5 min (blue)
Lysozyme reacted PEG for 30 min (red)



Changes in PEGylation with Reactant Ratios—Lysozyme




PEGylated Lysozyme—Deconvoluted MS Spectra



Expected Masses Correlate with Observed Masses

Protein	Expected Mass (kDa)	Observed Mass (kDa)
Lysozyme	14.3	14.3
Lysozyme + 1 PEG	16.6	16.6
Lysozyme + 2 PEG	18.9	18.9
Lysozyme + 3 PEG	21.2	21.2
Lysozyme + 4 PEG	23.5	23.5
Lysozyme + 5 PEG	25.8	25.8
Lysozyme + 6 PEG	28.1	Not Detected

Degree of PEGylation Tracked by Cation Exchange



Changes in PEGylation with Reactant Ratios—RNase A

Columns: ProPac SCX-10G Guard, 2 × 50 mm, ProPac SCX-10 Analytical, 2 × 250 mm

Mobile Phases: A) 20 mM MES, pH 6.1


B) 20 mM MES, 1M NaCl, pH 6.1

Gradient: 0–30% B from 0 to 30 min


Equilibration at 0% B for 10 min before injection

Flow Rate: 0.25 mL/min

Inj. Volume: 10 μ L



PEGylated RNase A – Deconvoluted MS Spectra



Expected Masses Correlate with Observed Masses

Protein	Expected Mass (kDa)	Observed Mass (kDa)
Ribonuclease A	13.7	13.4, 13.7
Ribonuclease A + 1 PEG	16.0	16.0
Ribonuclease A + 2 PEG	18.3	18.3
Ribonuclease A + 3 PEG	20.6	20.6
Ribonuclease A + 4 PEG	22.9	22.9
Ribonuclease A + 5 PEG	25.2	25.2
Ribonuclease A + 6 PEG	27.5	27.5
Ribonuclease A + 7 PEG	29.8	Not Detected

Conclusions

- PEGylated proteins were separated using SCX chromatography, revealing numerous isomers.
- These isomers were separated by the degree of PEGylation as confirmed by MS. PEGylation can be monitored as a function of reaction conditions, allowing reaction optimization for the desired product.
- SCX chromatography is shown to be a powerful method for characterizing modified protein mixtures with multiple positional isomers.

Thank you for your attention!