

Glycan Structural Elucidation On A Novel Quadrupole Dual Cell Linear Ion Trap Orbitrap Hybrid Mass Spectrometer

Julian Saba¹, Katie Southwick,² Shannon Eliuk,¹ Vlad Zabrouskov,¹ Sergei Snovida³

¹Thermo Fisher Scientific, San Jose, CA USA; ²Thermo Fisher Scientific, West Palm Beach, FL USA;

³Thermo Fisher Scientific, Rockford, IL, USA



Overview

Purpose: To demonstrate the use of HCD MS³ for glycan structural elucidation on a novel hybrid mass spectrometer, based on mass resolving quadrupole, Thermo Scientific™ Orbitrap™ analyzer, collision cell, linear ion trap architecture.

Methods: HCD MS² and MS³ spectra of permethylated chicken ovalbumin glycans were acquired on a Thermo Scientific™ Orbitrap Fusion™ Tribrid™ mass spectrometer. Structural elucidation was performed using SimGlycan® software

Results: The combination of permethylation, HCD MS³ and SimGlycan software enabled successful identification and differentiation of structural isomers of chicken-ovalbumin-released glycans

Introduction

Glycans in glycoproteins are involved in a wide-range of biological and physiological processes including recognition and regulatory functions, cellular communication, gene expression, cellular immunity, growth and development. Mass spectrometry (MS) has emerged as a powerful tool for glycan structural elucidation. The use of permethylation in combination with multistage fragmentation (MSⁿ) is critical to the success of this approach. Only analysis by MSⁿ truly characterizes glycans as it allows identification of heterogeneity, branching, linkages and resolution of isobaric structures which are otherwise indistinguishable in MS² spectra. Traditionally, MSⁿ has been restricted to low-energy, collisional-induced dissociation (CID) in linear ion trap mass spectrometers, thereby, requiring multiple stages of fragmentation (MS⁶, MS⁷...) for structural elucidation. Here we demonstrate for the first time the use of higher-energy collisional dissociation (HCD) MSⁿ for glycan structural elucidation on a novel instrument, the Orbitrap Fusion™ Tribrid™ mass spectrometer based on a mass resolving quadrupole, Orbitrap analyzer, collision cell, linear ion trap (Q-OT-qIT) architecture (Figure 1). The primary advantage of HCD fragmentation is the production of glycosidic, cross-ring, and internal double cleavage ions at the MS² level, where branching, linkage and resolution of isobaric structures are derived from the latter two types of ions. The availability of HCD MSⁿ enables comprehensive glycan structural elucidation at much lower MSⁿ stages.

Methods

Sample Preparation

Ovalbumin (1 mg, Sigma) was reduced, alkylated and digested overnight with trypsin in of 25 mM ammonium bicarbonate buffer (pH=8) at 37 °C. PNGase F solution (3 µL, Roche) was added to 200 µL of digested sample and the mixture was incubated for another 16 hours at 37 °C. The released glycans were separated from the peptides using a Sep-Pak® C18 cartridge (Waters). The Sep-Pak C18 was conditioned by washing with acetonitrile, followed by water. PNGaseF digested sample was loaded onto the cartridge and the released glycans were eluted with 1% ethanol while the peptides remained bound to the Sep-Pak C18. The released ovalbumin oligosaccharides were first purified using a porous graphite carbon column (PhyNexus) and then permethylated using an in-house protocol.

Mass Spectrometry


All MS experiments were performed using an Orbitrap Fusion Tribrid mass spectrometer via direct infusion into the nano-electrospray source. The mass spectrometer settings and SimGlycan software [1,2] version 4.50 (PREMIER Biosoft International) search parameters are listed in Tables 1 and 2.

Table 1. Mass Spectrometer Settings

Source	HESI	Isolation Width	3
Capillary Temperature	250 °C	Collision Energy	35
S-lens RF Level	60 %		
Source voltage [kV]	3.8		
Full MS Mass Range	700-1600 (<i>m/z</i>)		
MS Resolution	60000 @ <i>m/z</i> 200		
MS/MS Resolution	60000 @ <i>m/z</i> 200		

Table 2. SimGlycan Software 4.50 Settings

Ion Mode	Positive	Class	Glycoprotein
Adducts	Sodium	SubClass	N-Glycan (All)
Precursor m/z Error Tolerance	10 ppm	Biological Source	Chicken, Ovalbumin
Spectrum m/z Error Tolerance	0.01 Da	Pathway	Unknown
Chemical Derivatization	Permethylated	Search Structure	All
Reducing Terminal	Free	Glycan Type	All

FIGURE 1. Schematic representation of the novel Orbitrap Fusion Tribrid mass spectrometer architecture.

Results

HCD fragmentation provides much more informative spectra for glycan analysis compared to low-energy CID fragmentation available on linear ion trap mass spectrometers, (Figure 2). Targeting glycans with HCD fragmentation results in the production of cross-ring fragment ions and internal double-cleavage ions which, even with permethylation, can be lacking with low energy CID fragmentation. The ability to generate these types of ions at the MS² level provides sufficient information in a lot of cases to successfully elucidate glycan structures.

Though HCD can produce much more informative fragmentation at the MS² level, differentiating structural isomers can be an issue as fragments from mixed spectrum can complicate spectral assignment. MSⁿ would still be required for differentiation of these isomers. Currently the acquisition of HCD fragmentation is still limited to MS² on commercial mass spectrometers. However, we have recently built a new quadrupole, dual cell linear ion trap, Orbitrap hybrid mass spectrometer that has a novel architecture and ion transfer path that enables HCD MSⁿ analysis.

Structural elucidation using HCD MSⁿ was initially tested on glycans released from chicken ovalbumin. Since the glycan content of ovalbumin has been characterized in depth,^[3] it was an ideal system to examine the capabilities of the novel Orbitrap Fusion Tribrid mass spectrometer. Figure 3 shows the MS profile of permethylated glycans derived from ovalbumin acquired on this instrument. Data-dependent MS² spectra were acquired on all precursors with a charge state greater than two ($Z \geq 2$) and for each MS² spectra subsequent top 20 MS³ spectra were acquired.

The MSⁿ data interpretation workflow is as follows. The acquired MS² data is brought into SimGlycan software for automatic compositional identification. Based on the criteria selected (Table 2), SimGlycan software searches its database for glycan matches. For example the MS² spectrum (Figure 4) for the glycan at m/z 1046.511 (Figure 2) is brought into SimGlycan for identification. Based on the MS² fragmentation pattern SimGlycan interprets the spectrum as most likely being a hybrid glycan with a bisecting GlcNAc, as this is the glycan ranked as the top structure in SimGlycan search result (Figure 5).

It should be pointed out that the HCD MS² and MS³ spectra generated are information rich (Figure 6), containing informative glycosidic and cross-ring fragment ions. An additional advantage of HCD fragmentation is that the ions are measured within the Orbitrap mass analyzer with high-resolution, accurate-mass (HR/AM). This allows for differentiation of near mass fragment ions, which we observed to be very useful for correctly assigning branching and linkage ions.

Though we were able to confidentially identify the hybrid glycan with a bisecting GlcNAc (rank 1), there were additional structural isomers present. Closer examination of the MS² spectrum for the ion at *m/z* 1046.511 reveals numerous unidentified peaks that we could not attribute to the top ranked glycan. The acquisition of MS³ spectra on these peaks enabled identification of additional isomers. For example, the fragment ion observed at *m/z* 1606.789 in the MS² spectrum (Figure 4) indicates the loss of Gal-GlcNAc, which is most likely to occur from the non-reducing end of either asialyl digalactosyl biantennary glycan (ranked 4) or the other hybrid glycan with bisecting GlcNAc (ranked 3). The acquisition of MS³ spectrum for the ion at *m/z* 1606.789 provided additional fragment ions unique to the asialyl digalactosyl biantennary glycan that confirmed the presence of this structure at *m/z* 1046.511.

Overall, we were able to identify three structural isomers at *m/z* 1046.511. An additional hybrid glycan with the bisecting GlcNAc (ranked 3) was also identified. We had previously characterized these structural isomers on a Thermo Scientific™ Velos Pro™ mass spectrometer using CID MSⁿ.^[1] In order to differentiate the structural isomers, it required MS⁷, MS⁸ stages of fragmentation with CID MSⁿ.

Conclusion

- The novel Orbitrap Fusion Tribrid mass spectrometer based on mass resolving quadrupole, Orbitrap, collision cell, linear ion trap (Q-OT-qIT) architecture enables HCD MSⁿ.
- Additional glycan fragmentation pathways are accessible with this mass spectrometer as it has the ability to use any fragmentation mode, at any stage of MSⁿ analysis.
- The combination of permethylation, HCD MS³ and SimGlycan software enabled successful identification and differentiation of structural isomers of chicken ovalbumin released glycans.

References

1. Saba, J.; Apte, A.; Meitei, N.S.; Viner, R., Application Note 516: Automated Glycan Structural Isomer Differentiation Using SimGlycan Software.
2. Apte, A.; Meitei, N. S., Bioinformatics in glycomics: glycan characterization with mass spectrometric data using SimGlycan. *Methods Mol. Biol.* **2010**, 600, 269-281.
3. Harvey, D. J.; Wing, D. R.; Kuster, B.; Wilson, I. B., Composition of N-linked carbohydrates from ovalbumin and co-purified glycoproteins. *J Am Soc Mass Spectrom* **2000**, 11, (6), 564-71.
4. Domon, B.; Costello, C.E. A systematic nomenclature for carbohydrate fragmentations in FAB-MS-MS spectra of glycoconjugates. *Glycoconj J* **1988**, 5, 397-409.

www.thermoscientific.com

©2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. SimGlycan is a trademark of PREMIER Biosoft International. Sep-Pak is a trademark of Waters Corporation. All other trademarks are the property of Thermo Fisher Scientific, Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa-Other +27 11 570 1840
Australia +61 3 9757 4300
Austria +43 1 333 50 34 0
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China +86 10 8419 3588
Denmark +45 70 23 62 60

Europe-Other +43 1 333 50 34 0
Finland/Norway/Sweden
+46 8 556 468 00
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9434
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Russia/CIS +43 1 333 50 34 0
South Africa +27 11 570 1840



Thermo Fisher Scientific,
San Jose, CA USA
is ISO 9001:2008 Certified.

Spain +34 914 845 965
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC
Part of Thermo Fisher Scientific