

Identification and Quantification of 22 Common Anions in Pharmaceuticals in a Single Run Using HPIC with Suppressed Conductivity and Charge Detection

Hua Yang, Linda Lopez; Thermo Fisher Scientific, Sunnyvale, CA, USA



Overview

- An anion analysis method was developed to identify and quantify multiple counter and impurity ions
- The method described uses a capillary HPIC system with suppressed conductivity (CD) and charge (QD) detectors and a Thermo Scientific™ Dionex™ IonPac™ AS11HC-4µm capillary column. The method separates 22 common pharmaceutical anions in a single run.
- The counter ion profile of an allergy drug was analyzed using the method. Chloride was quantified with high sensitivity (LOQ = 0.004 mg/L) and a large linear range (0.1 to 500 mg/L).

Introduction

Ion analysis is important for the pharmaceutical industry because many active pharmaceutical ingredients (APIs) exist in their salt form. Pharmaceutical products are strictly regulated by the United States Food and Drug Administration (U.S. FDA) and other regulatory agencies, and must be tested for composition to verify their identity, strength, quality and purity.¹ The importance of ion analysis for the pharmaceutical industry has increased in recent years, not only because ion analysis is required for the final products in their salt form, but also ion analysis is becoming more critical for early stage drug development. Identification and quantification of ions is used in material quality control and counter ion selection during early stage drug development. The monitoring of material quality has become increasingly critical due to increased outsourcing of API synthesis. The APIs received from vendors may be contaminated with different counter ions during their synthesis, and therefore, have different physicochemical (solubility, crystallinity, hygroscopicity) and pharmacokinetic properties; Additionally, selecting a proper counter ion to improve API solubility and stability has become an increasing critical step in formulation development due to a trend that new drug-like compounds are less and less aqueous soluble.

Conventional HPLC is not suitable for these ion analysis, because they are typically small ions that lack a chromophore and are not retained by HPLC column. Ion chromatography (IC) with suppressed conductivity detection (CD) is an established and sensitive method used for ion determination.^{2,3} IC has been used previously in U.S. Pharmacopeia (USP) methods for single ion analysis of final pharmaceutical products. However ion analysis method is still needed to provide both counter and impurity ion profiles for APIs and formulations, as demanded by recent pharmaceutical developments.

In this application, 26 anions, which cover the majority counter and impurity ions in pharmaceutical samples, were selected for study with the objective to identify and quantify the analytes in a single run. The counteranions were separated on a high-capacity and high efficiency, 4 µm particle ion-exchange columns, facilitated on a Reagent-Free™ IC (RFIC™) high-pressure™ capable capillary IC system (HPIC™). The counteranions were detected by suppressed conductivity detection (CD) and the Thermo Scientific Dionex Charge Detector (QD). The resulting method is powerful and easy-to-use tool applicable for the pharmaceutical industry.

Experimental

Standards and Sample Preparation

Stock standards: Twenty-six individual stock standards (~1000 mg/L) were prepared by dissolving reagents for: chloride, gluconate, acetate, glycolate, formate, pyruvate, glucuronate, nitrate, bromide, glutarate, succinate, malate, tartrate, malonate, benzoate, maleate, sulfate, fumarate, phosphate, citrate, tosylate, benzenesulfonate, lactate, fluoride, nitrite, and trifluoroacetate in 18 MΩ-cm resistivity deionized water.

Mixed standard: Prepare the mixed standard from the 1000 mg/L anion stock standards and deionized water.

Chloride calibration standards: Dilute the 1000 mg/L chloride stock solution appropriately with deionized water to make nine chloride calibration standards ranging from 0.1 to 500 mg/L.

Samples: An over-the-counter (OTC) allergy drug with API as a chloride salt was purchased from a local pharmacy. Samples were prepared by dissolving one tablet in 1000 mL of deionized water and then filtering with a 0.20 µm IC syringe filter prior to injection to remove insoluble particles. Some filtered samples were diluted with deionized water.

Ion Chromatography Conditions

Instrument:	Thermo Scientific Dionex ICS-5000+ HPIC capillary system
Column:	Dionex IonPac AS11-HC-4µm, 0.4 × 250 mm
Eluent Source:	Thermo Scientific Dionex EGC-KOH capillary cartridge
Gradient:	1.5–2 mM KOH (0–5 min), 2–8 mM (5–8 min), 8–16 mM (8–26 min), 16–70 mM (26–32 min), 70 mM (32–38 min).
Flow Rate:	0.0150 mL/min
Inj. Volume:	0.40 µL
Column Temp.:	30 °C
IC Cube Temp.:	15 °C
Detection:	CD: Suppressed conductivity detection, AutoSuppression™ Thermo Scientific™ Dionex™ ACES™ 300 Anion Electrolytic Suppressor, recycle mode QD: Charge detection, 6V

Data Analysis

Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System (CDS) software.

Figure 1 shows the schematic of the flow diagram of the HPIC system used for this application. It is a high-pressure capillary IC system (HPIC) with CD / QD detectors. The capillary IC components are configured in close proximity in Thermo Scientific™ Dionex™ IC Cube™ module (Figure 2) to minimize the flow path.

FIGURE 1. Schematic of capillary HPIC system with dual detectors.





FIGURE 2. Thermo Scientific Dionex IC Cube module.




Results and Discussion

Twenty-six of most commonly used counter ions and commonly seen impurity ions in pharmaceutical samples—were included in this study. The Dionex IonPac AS11-HC-4 μ m capillary column was selected for its ability to separate various organic and inorganic acids in a wide range of products including-fruit juices and wine.⁴ The 4 μ m resin particle technology has further improved the resolution and increased the column efficiency. As can be seen in Figure 3, 22 of 26 organic and inorganic anions were separated in a single run in less than 40 minutes. Four anions coelute with the other ions of interest, and therefore removed from the mixed standard.

The QD detector is an excellent orthogonal complement for the CD detector. Because each ion had a characteristic retention time, CD response, and QD response (Figure 3), all ions were identified and quantified with confidence using the two detectors.

FIGURE 3. Chromatogram of an anion standard mix.



Chloride is the most prominent counter ion used in pharmaceuticals. The application example reported here uses the method to analyze an allergy drug in which both APIs are chloride salts.

Figure 4-1 shows the chromatograms of a water blank, an allergy drug tablet dissolved in 1000 mL water, and its 5-fold dilution. The results confirm that the drug has a clean counter ion profile with mainly chloride and trace amount of acetate, nitrite, nitrate and sulfate.

Figure 4-2 compares the CD and QD chromatograms of the undiluted sample. An additional unknown anion (peak 3) detected by QD illustrates the advantages of having two detectors that are based on different technologies. For some analytes, CD can be used to detect and quantify them, but with QD they are opaque or minimally detected and vice versa.

Figure 4-1 Chromatograms using suppressed conductivity (CD) detection.



Conclusion

- A versatile and easy to use anion analysis method was developed for the pharmaceutical industry which is suitable for the study of both counter and impurity ion profiles for APIs and formulations.
- Using a capillary HPIC system with CD and QD detectors, and a Dionex IonPac AS11HC-4µm capillary column:
 - Twenty-two common pharmaceutical anions were separated in a single run.
 - Multiple counter ions in drug products were easily identified and quantified with confidence, which facilitated the profiling of this allergy drug.
 - For chloride, this method had a low limit of quantification (0.004 mg/L) and a large linear range (0.1 to 500 mg/L).

References

1. FDA, CFR Title 21, Foods and Drugs, Chapter 1, Part 211 *Current Good Manufacturing Practices for Finished Pharmaceuticals*
<http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfcfr/CFRSearch.cfm?fr=211.110> (accessed April 15, 2014).
2. Dionex (now part of Thermo Scientific) Application Note 116: Quantification of Anions in Pharmaceuticals. Sunnyvale, CA, 2002. [Online]
http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/Application%20&%20Technical%20Notes/Chromatography/GC%20HPLC%20and%20UHPLC%20Columns%20and%20Accessories/Chromatography%20Column%20Accessories/4079-AN116_LPN0924-01.pdf (accessed April 15, 2014).
3. Dionex (now part of Thermo Scientific) Application Note 164: Assay for Citrate and Phosphate in Pharmaceutical Formulation Using Ion Chromatography. Sunnyvale, CA, 2004. [Online]
http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/Application%20&%20Technical%20Notes/Chromatography/GC%20HPLC%20and%20UHPLC%20Columns%20and%20Accessories/Chromatography%20Column%20Accessories/49183-AN164_LPN1643.pdf (accessed April 14, 2014).
4. Dionex (now part of Thermo Scientific) Application Note 273: Higher Resolution Separation of Organic Acids and Common Inorganic Anions in Wine. Sunnyvale, CA, 2011. [Online]
<http://www.thermoscientific.com/content/dam/tfs/ATG/CMD/CMD%20Documents/Application%20&%20Technical%20Notes/Chromatography/Liquid%20Chromatography/Liquid%20Chromatography%20Accessories/110622-AN273-IC-OrganicAcids-Anions-Wine-15Jun2011-LPN2727-01.pdf> (accessed April 15, 2014).

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, Sunnyvale, CA
USA is ISO 9001:2008 Certified.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand