

Nitrogen, Carbon and Sulfur Determination in Paper by Flash Combustion

Liliana Krotz and Guido Giuzzi
Thermo Fisher Scientific, Milan, Italy



Overview

Purpose: To show the characterization of paper samples by Organic Elemental Analysis (OEA).

Methods: Paper samples were analyzed using an elemental analyzer with an automatic autosampler.

Results: Data collected of nitrogen, carbon and sulfur from different paper samples are discussed to assess the performance of the OEA analyzer.

Introduction

In the production process of paper, elemental composition is periodically monitored and tested for the characterization of raw and final products. Nitrogen and carbon are the most important parameters in quality control whilst the sulfur content is an indication of impurities present in the materials.

The Thermo Scientific™ FLASH 2000 analyzer (Figure 1) permits the fast, quantitative determination of elements in paper materials without any sample pre-treatment. The system, which is based on the dynamic combustion of the sample, provides automatic and simultaneous nitrogen, carbon and sulfur determination in a single analysis run.

FIGURE 1. FLASH 2000 Elemental Analyzer



Method

The sample is weighed in a tin capsule and introduced into the combustion reactor via the Thermo Scientific™ MAST™ 200R autosampler together with the proper amount of pure oxygen. In the NCS configuration, the gases produced after combustion of the sample are carried by a helium flow through a copper layer, and then swept through a water trap and a GC column that separates the combustion gases which are finally detected by a Thermal Conductivity Detector. Total run time is 10 minutes (Figure 2).


For nitrogen only determination, the gases produced after combustion are carried by a helium flow to a second reactor filled with copper, then swept through CO₂ and H₂O traps, then through a GC column and finally sensed by a Thermal Conductivity Detector. Total run time is less than 5 minutes (Figure 3).

A complete report is automatically generated by the Thermo Scientific™ Eager Xperience dedicated data handling software system and displayed at the end of the analysis.

Analytical Conditions – NCS determination

| | |
|------------------------|---------------|
| Reactor Temperature: | 950 °C |
| Oven Temperature: | 65 °C |
| Helium Carrier Flow: | 140 ml/min |
| Helium Reference Flow: | 100 ml/min |
| Oxygen Flow: | 250 ml/min |
| Oxygen Injection End: | 5 sec |
| Sampling Delay Time: | 12 sec |
| Run Time: | 600 sec |
| Standard: | 2 - 3 mg BBOT |
| Sample Weight: | 3 - 4 mg |


FIGURE 2. NCS configuration.



Analytical Conditions – NCS determination

| | |
|---------------------------------|-------------------------------------|
| Combustion Reactor Temperature: | 950 °C |
| Reduction Reactor Temperature: | 840 °C |
| Oven Temperature: | 50 °C |
| Helium Carrier Flow: | 140 ml/min |
| Helium Reference Flow: | 100 ml/min |
| Oxygen Flow: | 300 ml/min |
| Oxygen Injection End: | 20 sec |
| Sampling Delay Time: | 10 sec |
| Run Time: | 360 sec |
| Standard: | 40 – 50 mg Aspartic Acid (10.52 %N) |
| Sample weight: | 90 – 100 mg |

FIGURE 3. Nitrogen determination



Results

Different paper samples, cut into small pieces were chosen to assess the system. Samples were analyzed several times to evaluate the reproducibility of the method.

Table 1 shows the NCS data of paper samples. Whilst in the NCS configuration, the system was calibrated with 3 – 4 mg of BBOT* standard using K factor as calibration method. The sample weight of paper was 3 - 4 mg. Samples were analyzed with the addition of Vanadium Pentoxide for a complete conversion of sulfur. Figure 4 shows a typical NCS chromatogram.

*BBOT: 2,5 bis(5-tert-butyl-benzoxazol-2-yl)thiophene: 6.51 %N, 72.53 %C, 7.44 %S

TABLE 1. NCS determination.

| Sample | N % | RSD % | C % | RSD % | S % | RSD % |
|---------|--------|--------|---------|--------|--------|--------|
| Paper A | 1.8279 | 1.2421 | 42.0787 | 0.1086 | 0.0437 | 0.9556 |
| | 1.8310 | | 42.1599 | | 0.0439 | |
| | 1.8690 | | 42.1559 | | 0.0431 | |
| Paper B | 0.0282 | 3.0270 | 42.3373 | 0.0467 | 0.0394 | 4.3372 |
| | 0.0298 | | 42.3655 | | 0.0422 | |
| | 0.0284 | | 42.3274 | | 0.0390 | |
| Paper C | 0.9866 | 1.1261 | 42.4505 | 0.0517 | 0.0422 | 0.9592 |
| | 1.0042 | | 42.4331 | | 0.0425 | |
| | 0.9835 | | 42.4767 | | 0.0417 | |
| Paper D | 1.2053 | 1.0200 | 42.3358 | 0.0600 | 0.0360 | 0.5587 |
| | 1.2293 | | 42.3647 | | 0.0356 | |
| | 1.2119 | | 42.3141 | | 0.0358 | |
| Paper E | 1.2529 | 0.1195 | 42.1291 | 0.2729 | 0.0374 | 1.5851 |
| | 1.2521 | | 42.2987 | | 0.0363 | |
| | 1.2550 | | 42.3493 | | 0.0372 | |

FIGURE 4. Typical NCS chromatogram




Table 2 shows the data of nitrogen determination of the electrostatic paper E in order to compare the results obtained previously by the NCS configuration. The instrument was calibrated with Aspartic acid (10.52 %N) as standard using K factor as calibration method. Table 3 shows the relative comparison, demonstrating that there are no significant differences in nitrogen percentages using different weight ranges. This confirms the complete conversion of the element during combustion.

TABLE 2. Nitrogen data of Electrostatic Paper E.

| Weight (mg) | N % | Average N % | RSD % |
|-------------|------|-------------|-------|
| 99.8 | 1.25 | 1.25 | 1.04 |
| 107.2 | 1.23 | | |
| 106.0 | 1.26 | | |
| 102.1 | 1.24 | | |
| 107.4 | 1.26 | | |

TABLE 3. Nitrogen comparison data of Electrostatic Paper E by NCS and N determination.

| NCS Configuration | | | | Nitrogen configuration | | | |
|-------------------|--------|-------------|--------|------------------------|------|-------------|-------|
| Weight (mg) | N % | Average N % | RSD % | Weight (mg) | N % | Average N % | RSD % |
| | | | | 99.80 | 1.25 | | |
| 3.958 | 1.2329 | | | 107.20 | 1.23 | | |
| 4.025 | 1.2521 | 1.2533 | 0.1195 | 106.00 | 1.26 | 1.25 | 1.04 |
| 3.969 | 1.2550 | | | 102.10 | 1.24 | | |
| | | | | 107.40 | 1.26 | | |

Figure 5 shows the curve calibration obtained for nitrogen only determination.

FIGURE 5. Nitrogen Curve Calibration.




Table 4 shows the data of nitrogen determination of other paper samples. Whilst in the nitrogen configuration, the system was calibrated with 50 - 100 mg of Aspartic Acid (10.52 %N) and the paper sample weight was 70 - 100 mg.


TABLE 4. Nitrogen Determination.

| Sample | N % | RSD % |
|---------|----------------------------|--------|
| Paper 1 | 0.3282 0.3368 0.3326 | 1.2932 |
| Paper 2 | 0.3334 0.3372 0.3344 | 0.5880 |
| Paper 3 | 0.2513 0.2501 0.2469 | 0.9119 |
| Paper 4 | 0.2959 0.3003 0.2963 | 0.8179 |
| Paper 5 | 0.2880 0.2930 0.2888 | 0.9263 |



Figure 6 shows a typical chromatogram of Nitrogen only determination

FIGURE 6. Typical Nitrogen chromatogram.



Conclusion

All data were obtained with a good reproducibility and no matrix effect was observed when changing the sample.

We demonstrate that the advantage of the FLASH 2000 analyzer lies in its ability to perform NCS determination in a single run, then, by changing the configuration and increasing the sample weight, it is possible to perform nitrogen only determination.

Using the elemental analyzer, it is also possible to characterize the different industrial paper applications according to the chemical concentration of the elements. This makes it possible to choose the most suitable paper recycling system, and perform a Life Cycle Assessment (LCA) on the environmental impact of the products from their production to their disposal.

This poster demonstrates that the FLASH 2000 OEA copes with all the demanding requirements of modern laboratories such as flexibility, accuracy, reproducibility, sensitivity, and automation.

www.thermoscientific.com

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific, Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

| | | | |
|--|---------------------------------------|--------------------------------------|-------------------------------------|
| Africa +43 1 333 50 34 0 | Denmark +45 70 23 62 60 | Japan +81 45 453 9100 | Singapore +65 6289 1190 |
| Australia +61 3 9757 4300 | Europe-Other +43 1 333 50 34 0 | Latin America +1 561 688 8700 | Spain +34 914 845 965 |
| Austria +43 810 282 206 | Finland +358 9 3291 0200 | Middle East +43 1 333 50 34 0 | Sweden +46 8 556 468 00 |
| Belgium +32 53 73 42 41 | France +33 1 60 92 48 00 | Netherlands +31 76 579 55 55 | Switzerland +41 61 716 77 00 |
| Canada +1 800 530 8447 | Germany +49 6103 408 1014 | New Zealand +64 9 980 6700 | UK +44 1442 233555 |
| China 800 810 5118 (free call domestic) 400 650 5118 | India +91 22 6742 9494 | Norway +46 8 556 468 00 | USA +1 800 532 4752 |
| | Italy +39 02 950 591 | Russia/CIS +43 1 333 50 34 0 | |



Thermo Fisher Scientific,
San Jose, CA USA
is ISO 9001:2008 Certified.

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific