

Achieve Superior Sample Integrity and Seal in a Bargain!

Three Great offers on Capping and Sealing Products



Seal in a Bargain!

Great offers on Capping and Sealing Products

Choose from a range of capping and sealing options for long-term storage, colorimetry, fluorescence, PCR and more. Ensure sample security while speeding sample processing... and get a great price!

- Purchase a Capit-All[®] or Capit-All IS and Thermo Scientific Matrix ScrewTop Tubes to receive a **25% discount on your order!**
- Purchase a Thermo Scientific ALPS 3000 Heat Sealer and receive a **33% discount and a roll of heat seals FREE!**
- Purchase a Thermo Scientific ALPS 50V Microplate Heat Sealer and receive a **Thermo Scientific ALPS 25 Sealer FREE!**

Purchase a Capit-All® or Capit-All IS and Thermo Scientific Matrix ScrewTop Tubes to receive a 25% discount on your order!

Save on capping in high-throughput environments

Cap or decap a full rack of 96 tubes in approximately 10 seconds!

Increase throughput with precision, ease and peace of mind.

Compact, intuitive Capit-All decappers are designed specifically to work with Matrix® ScrewTop tubes. Seal integrity is maintained by precise torque control.

Ideal for the highest throughput environments, the Capit-All IS instrument easily can be integrated via RS32.

Choose either instrument:

Cat. No.	Description
4111MAT	Capit-All Capper/Decapper
411MAT-IS	Capit-All IS Automation Friendly Capper/Decapper

And any one or more of the following items:

Cat. No.	Description	Quantity
3741*	Matrix 1.0 ml ScrewTop Tubes	5 Latch Racks of 96 tubes/case
3742	Matrix 1.0 ml ScrewTop Tubes w/o caps	5 Latch Racks of 96 tubes/case
3744*	Matrix 500 µl ScrewTop Tubes	5 Latch Racks of 96 tubes/case
3745	Matrix 500 µl ScrewTop Tubes w/o caps	5 Latch Racks of 96 tubes/case
4477*	Matrix ScrewTop Cap Tray	5 Trays/case
4906	Matrix ScrewTop Cap Tray, empty	5 Trays/case

All ScrewTop Tubes are Sterile

*Options with colored caps are available for these configurations. Contact your sales representative for details or visit www.thermoscientific.com/matrix.



Purchase a Thermo Scientific ALPS 3000 Heat Sealer and receive a **33% discount** and a roll of heat seals **FREE!**

Save on microplate heat sealing in high-throughput applications

Seal up to 10 plates per minute!

The easy-to-use, flexible ALPS™ 3000 features touch screen operation for benchtop sealing. The unit integrates with robotics without special software. The automatic system shut - off eliminates product waste, system hold-ups and potential injury.

Seal options are available for PCR, colorimetry, fluorescence, long-term storage, piercing and resealing and more.



Purchase an ALPS 3000 Sealer with a 33% discount:

Cat. No.	Description
AB-3000	ALPS 3000

And select your free roll of seals:

Cat. No.	Description
AB-3559	Thermo Seal
AB-3686	Clear Seal Strong
AB-3720	Easy Pierce 20 µm
AB-3738	Easy Pierce
AB-3739	Easy Peel
AB-3797	Clear Seal
AB-3799	Diamond Seal



Purchase a Thermo Scientific ALPS 50V Microplate Heat Sealer and receive a Thermo Scientific ALPS 25 Sealer FREE!

Get two for the price of one! Save on microplate heat sealing for low-throughput applications

The ALPS 50V microplate heat sealer offers variable time and temperature settings for highly reproducible sealing and superior sample storage. Whether your application is PCR, sample archiving or compound storage, your research is protected.

Both sealers are compact, portable, and can be used with any plate or block height, including PCR plates. A wide variety of seals are available to suit most any plate storage need.


Purchase an ALPS 50V Sealer:

Cat. No.	Description
AB-1443	ALPS 50V

And receive your ALPS 25 Free:

Cat. No.	Description
AB-0384/110	ALPS 25





Act Now

to take advantage of these special offers!

1. To place your order and receive your discount, please contact your sale representative or call Customer Service at 800-345-0206.
2. To receive your free product once you've completed your purchase, visit: www.thermoscientific.com/redeemsamplestorage

You will be directed to a redemption form. Proof of purchase required. Free product shipped separately from order; allow 4 – 6 weeks for delivery. All offers valid January 18, 2010 to June 30, 2010. Orders must be placed directly with Thermo Fisher Scientific. One redemption for each submission made. Redemptions must be completed by August 15, 2010. Valid in U.S. and Canada only. Void where prohibited by law or company policy. Prices and specifications of products subject to change at any time. Not to be combined with any other offers.

© 2010 Thermo Fisher Scientific Inc. All rights reserved. Capit-All is a registered trademark of The Automation Partnership (TAP). All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
Europe: Tel: +44 (0) 161 496 2110 | email: matrix.eu.info@thermofisher.com
Asia: China: Tel: +86-21-6865 4588, Japan: Tel: +81-45-453-9220, India: Tel: +91-22-6716-2200,
Hong Kong & Other Asian Countries: Tel: +852-2885-4613, Australia: Tel: 1-300-735-292,
New Zealand: Tel: 0800-933-966 | email: matrix.ap.info@thermofisher.com
To order online: www.thermoscientific.com/Matrix

Thermo
SCIENTIFIC