


CUSTOM SEQUENCING COMB



End User _____

Phone _____

Company _____

Fax _____

Email _____

Owl Customer Service Rep _____

Proposal Number FOR OFFICE USE ONLY

Product: ☐ Owl Brand ☐ Other Brand

Multi channel pipette format? ☐ Yes ☐ No ☐ 1X or ☐ 2X

Quantity _____

OWL PRODUCT COMB LENGTH

Model No. L (cm)

ADJ1 13.3

ADJ2 17.1

ADJ3 30.48

S1S 17.1

S2S 30.48

S3S 30.48

S4S 17.1

SIMPLE - for use with Owl products only

Owl Unit/Model Number _____

Number of Teeth _____ Thickness _____ cm
(thickness of 0.08cm is .03 in / thickness of 0.015cm is 0.5in.)

COMPLEX - for use with Owl and other brands

Material: ☐ Teflon ☐ Polyester (use if thickness if less than .030)

Length _____ cm Height _____ cm

Number of Teeth _____ Tooth Height _____ cm Tooth Width _____ cm (0.3 min.)

Gap _____ cm (0.254 min.) Overhang _____ cm

OR Center to Center _____ cm

Slots: Si _____ cm Sd _____ cm Sw _____ cm

(Slots will be centered on upper area of comb)

$$L = [(Number\ of\ Teeth) \times (Tooth\ W)] + [Gap\ (Number\ of\ Teeth - 1)] + 2\ (Overhang)$$

PLEASE: Only fill out those dimensions that are different from a standard comb. We will fill in the rest.
THERE IS A MINIMUM CHARGE OF \$250.00 FOR ALL POLYESTER SHARKTOOTH SEQUENCING COMBS.