


# Technical Specifications

## Solenoid Valve, D2, MP-6A & ACS



EXTRA PARTS NEEDED:  
1) 440376 - Automatic Still Adapter Kit

NOTE:  
Cartridges are not included.  
Proper cartridge selection is based on inlet water quality and product water applications.

- NOTES:**
1. **CAUTION:** Refer to appropriate Barnstead Equipment manuals before installing, operating or servicing this equipment.
  2. **CAUTION:** Water from this unit is for a specific use. For best results, carefully follow the instructions in the appropriate Barnstead manuals.
  3. Installer to attach unit to wall with suitable fasteners (customer supplied).
  4. Do not mount electrical equipment or outlets below this unit.
  5. For additional information, see owners manual.

<b>Drawing No.</b>	SDST 21-93 (1741-92)
<b>Description</b>	Solenoid Valve, D2, MP-6A, ACS
<b>Scale</b>	Not to Scale
<b>Date</b>	02/26/04

Barnstead International • 2555 Kerper Blvd., P.O. Box 797 • Dubuque, Iowa 52004-0797 • 1-800-553-0039