



- Legend:
- 1 Barnstead GenPure Pro Ultrapure Water Purification System
 - 2 Holes are for optional wall mounting of water system with supplied hardware

Notes:

- 1 Dual dimensions are in metric and inch

ELECTRICAL REQUIREMENTS	100-240 V, 50/60 Hz, 2-1A, up to 5 ft from unit			
WATER CONNECTIONS	3/4" NPT with manual shut off valve			
MIN/MAX INLET PRESSURE	2 - 87 PSI (0.1 - 6 bar)			
RECOMMENDED FEED TEMPERATURE	2 - 4.0°C			
RECOMMENDED FEED WATER TYPE	Pure water (Type 2 water)			
DRAIN	An atmospheric drain must be available within 5 feet of the final mounting location			
Client / Town:	Project:			
Index	Amendment	Date	Name	Scale
a	gesdient	16.11.2014	M.A	
		Drawn	08.10.2014	Avadjeva
		Checked	08.10.2014	ODP
Thermo Electron LED GmbH Steinweg 3, D-55672 Niederdorf Tel. +49 (0)2021 106 99-0				
Dimension sheet GenPure Pro				
MB124				Sheet
				Sheets