


Connections from tank



**Legend:**

- 1 Double Pretreatment for Barnstead Pacific TII/RO water systems (Item No. 09.4.000)
- 2 Barnstead Pacific TII/RO pure water system
- 3 Barnstead 30L Reservoir for Pacific TII/RO water systems
- 4 OPTIONAL Hand Dispenser can be utilized up to 3m (108") from water system
- 5 Holes are for optional wall mounting of water system with supplied hardware
- 6 OPTIONAL Reservoir Wall Mounting Bracket (Item No. 06.5016)

**Notes:**

- 1 Dual dimensions are in metric and inch

ELECTRICAL REQUIREMENTS		100-240 V, 50/60 Hz, 2-1A, up to 5 ft from unit	
WATER CONNECTIONS		3/4" NPT with manual shut off valve recommended	
MIN/MAX INLET PRESSURE		15 - 85 PSI (0.1 - 6 bar)	
RECOMMENDED FEED TEMPERATURE		2 - 35°C	
RECOMMENDED FEED WATER TYPE		Tap/Potable Water	
DRAIN		An atmospheric drain must be available within 5 feet of the final mounting location	
Client / Town:	Project:	Date:	Name:
		07.10.2014	Avadjeva
		07.10.2014	ODP

<p><b>Thermo</b>                  SCIENTIFIC</p> <p>Thermo Electron LED GmbH                  Schwelm 3, D55672 Niedewald                  Tel. +49 (0)202 1 06 99-0</p>		<p>Dimension sheet</p> <p>Pretreatment, Pacific, 30L Tank</p> <p>MB118</p>	
Sheet	Sheet	Date	Name
		07.10.2014	Avadjeva
		07.10.2014	ODP
Scale			