


Connections from tank



| | |
|------------------------------|---|
| ELECTRICAL REQUIREMENTS | 100-240 V, 50/60 Hz, 2-1A, up to 5 ft from unit |
| WATER CONNECTIONS | 3/4" NPT with manual shut off valve recommended |
| MIN/MAX INLET PRESSURE | 15 - 85 PSI (0.1 - 6 bar) |
| RECOMMENDED FEED TEMPERATURE | 2 - 35°C |
| RECOMMENDED FEED WATER TYPE | Tap/Potable Water |
| DRAIN | An atmospheric drain must be available within 5 feet of the final mounting location |

| | | |
|----------------|----------------------|--------|
| Client / Town: | Project: | Scale: |
| Date: | Name: | |
| Drawn: | 11.09.2014, Avadjeva | |
| Checked: | 11.09.2014, OPP | |

| | | | |
|---|---------------|------------|-----|
| c | Vorbehandlung | 26.06.2015 | M.A |
| b | geändert | 26.11.2014 | M.A |
| a | geändert | 14.11.2014 | M.A |

| | | | |
|--|--|--|-------|
| | | Thermo Electron LED GmbH Schwelm 3, D-55672 Niederdorf Tel. +49 (0)202 1 06 99-0 | |
| Dimension sheet Pretreatment, Pacific, 60L Tank | | Sheet | MB117 |
| | | Sheets | |

- Legend:
- 1 Double Pretreatment for Barnstead Pacific TII/RO ater systems (Item No. 09.4.000)
 - 2 Barnstead Pacific TII/RO pure water system
 - 3 Barnstead 60L Reservoir for Pacific TII/RO water systems
 - 4 OPTIONAL Hand Dispenser can be utilized up to 3m (108") from water system
 - 5 Holes are for optional wall mounting of water system with supplied hardware
 - 6 OPTIONAL Reservoir Wall Mounting Bracket (Item No. 06.5016)
- Notes:
- 1 Dual dimensions are in metric and inch